

# GLOWSKIN O<sup>+</sup>

VER:1.0

## USER MANUAL



---

## Warning


Dear user:

Thanks for using our company's products, in order to fully make use of the products, we sincerely suggests follow:


- ◆◆ Please read the instruction in detail and keep it well for further reading and consulting.
- ◆◆ Please follow the instruction to fit and operate the machine correctly.
- ◆◆ Please do not remove and change any accessories of the machine.
- ◆◆ It forbids others to open and remove the mainframe except the authorized people in the company.
- ◆◆ It forbids others to change and maintain the machine except the professional people.
- ◆◆ Please contact us if there is any trouble during operation.

**This product does not belong to medical,  
please don't use for disease treatment!**

## Parts Introduction



1. Induction card location
2. Touch screen
3. Vacuum adjust valve
4. Water flowing adjust valve
5. Cool & Heat work handle jack1
6. Cool & Heat work handle jack2
7. BOHR effect C&O handle jack
8. 4D-Circle RV handle jack
9. Sprayer jack
10. Aqua peeling handle water hose jack
11. Aqua peeling handle vacuum hose jack
12. Big Micro vibration work handle jack
13. Small Micro vibration work handle jack
14. Aqua peeling handle filter
15. Accessories frame
16. S4 liquid bottle
17. S3 liquid bottle
18. S2 liquid bottle
19. S1 liquid bottle
20. Dust filter
21. General power switch
22. Power supply line socket
23. Water recycling bottle
24. Water recycling bottle outlet



Any part of instrument made of PMMA material (for example: accessories shelf, panel, transparent handles, etc.). Do not use alcohol or corrosive solvent to wipe clean.

---

---

## Accessories



BOHR effect C&O handle



4D-Circle RV handle



Sprayer



Cool & Heat work handle



Small Micro vibration work handle



Big Micro vibration work handle



Aqua peeling clean handle



Aqua peeling clean cap

## Basal Operation

1.Put accessories well and connect them to the host.

(1)Insert the plug of BOHR effect C&O handle into[7] jack, tighten the nut and put it into frame.

(2)Insert the plug of 4D-Circle RV handle into [8] jack, tighten the nut and put it into frame.

(3)Insert the plug of sprayer into [9] jack and it into accessories frame.

(4)Insert the plug of Aqua peeling handle thin hose(water) into [10] jack.

(5)Insert the plug of Aqua peeling handle thick hose(vacuum) into [11] jack.

(6)Insert the plug of Cool & Heat work handle into [5][6] jack.

(7)Insert the big Micro vibration work head into [12] jack.

(8)Insert the small Micro vibration work head into [13] jack.

(9)Pour distilled water into the S4 cleaning bottle and install it into position [15].

(10)Pour liquid into the S3 cleaning bottle and install it into position [16].

(11)Pour liquid water into the S2 cleaning bottle and install it into position [17].

(12)Pour liquid water into the S1 cleaning bottle and install it into position [18].


(13)Install the Water recycling bottle in position [22].

2.Insert the plug of Power supply line into [21] socket, turn on Power switch [20]. Input system default password: "000000" , the system enter into (Fig 1) shows.

3.In the interface as (Fig 1):

(1)Click on "Basic setting" to select language, set date and time. **Attention: The "Advanced setting" is only for manufacturer, Users do not use.**

(2)Click on "item" to enter operating interface.



(Fig 1)

## Aqua Peeling

1. Click on “Aqua peeling” in (Fig 1) to enter interface as (Fig 2).

2. Clean the skin.

3. Pour application of liquid products into the liquid bottle.

S1 AHA Dermabrasion skin, activate collagen, strengthen skin nutrition and transport, Oxidation resistance.

S2 BHA acne, blackhead, sebum to grain, deep cleaning, sterilization, tighten pores.

S3 HYDRO vitamin B supplement, EGF and so on other nutrients.

S4 Pour distilled water.

3. Set work time. The system default work time is 30 minute. Press “+” and “-” of “time” to set. The regulation range is 1 to 120 minutes.

4. According to treatment part to choose big or small Aqua peeling clean cap, connect cap and Aqua peeling clean handle together.

5. Adjust flux. Turn water flowing adjust valve[3] knob to adjust Aqua peeling clean handle flux, turning it in clockwise direction to become large, turning it in count clockwise direction to become small.

6. Adjust vacuum intensity. Press the Aqua peeling clean cap with your finger before adjust, and then turn [2] Vacuum adjust valve knob the pressure is strong while turning it to the clockwise. turning it in count clockwise to become weak. The beautician should consider the nursing part and make sure the guests can accept the pressure.

7. Click on “Start”, the system begin working. Press “Pause” to stop working.

8. Beautician hold the Aqua peeling clean handle on skin for a while until the water flow into the cap, and then drop some water on skin by points, to make the skin surface wet, and then apply cleansing product on the skin, and move the it on skin back and forth.

9. The system will pause automatically if setting time is over. press any key to continue.

10. After finished operating, click on “S4”, Put the silicone sheet on the cap, adjust the vacuum and water flow to biggest, start the aqua peeling, the system begin clean hose.

11. Click on “Return” to back to (Fig 1).


### Cautions:

1. Do not use on wound skin.

2. Pour out water in waste recycling bottle in time to avoid damage caused by overflowing.

3. Different size of cap can be flush with water directly, then soaked with 75% alcohol.

**4. Turn off the machine. Clean water recycling bottle everyday. adjust Vacuum valve to min, the handle do not block. insert the connector to outlet, the water flowing.**



(Fig 2)



## BOHR Effect C&O

1. First, swipe the induction card under carbon oxygen head package to check "BOHR Effect C&O" function. click on "BOHR Effect C&O" in (Fig 1) to enter interface as (Fig 3).

2. The system default function is "Revive", Click on "Bright" and "Massage" to select other function.

3. Set LED color of handle. Click on "RED" or "Blue". The indicator on handle will light up. If choose "RED" and "Blue", the indicator color is purple.

4. Set intensity. The system default intensity is 10%. Click on "-" and "+" to choose. "Revive" adjust range is 10% to 30%. "Bright" adjust range is 10% to 60%. "massage" adjust range is 10% to 100%.

5. Set work time.

6. Connect Carbon oxygen head on BOHR Effect C&O handle.

7. Apply gels on operating parts. Click on "Start", The system in pre-start states. If you need to pause, Click on "Pause".

8. Beautician press "start" key on handle, The beautician holds the handle and moves slowly in circle. you can press "+" and "-" in handle to adjust intensity. The function is the same as interface.

9. The system will pause automatically if setting time is over. press any key to continue.

10. Click on "Return" to back to (Fig 1).

### Cautions:

1. Do not use on wound skin.

2. Carbon oxygen heads are consumables, after used, please replace in time.



(Fig 3)

## 4D Circle RV

1. Click on “4D circle RV” in (Fig 1) to enter interface as (Fig 4).
2. The system default operating parts is “Face”. Click on “Neck”, “Decollete”, “Hands” to choose deferent parts.
3. Set Power level. The system default power level is “Low”, Click on “Medium” and “High” to choose deferent level .



(Fig 4)

4. Set photon mode. The system default mode is P0(Off). Click on “-” and “+” of “photon mode” to set P1 is continuous light. P2 to P4 is pulse light.
5. Set intensity. The system default work intensity is 10%. Click on “-” and “+” of “intensity” to regulate. Low level is 10% to 30%. Medium level is 10% to 60%. High level is 10% to 100%
6. Set work time.
7. Apply gels on operating parts. Click on “Start” , The system in pre-start states. If you need to pause, Click on “Pause”.
8. Beautician press “start” key on handle, The beautician holds the handle and moves slowly in circle. Do not stop at the same parts without operating. you can press “+” and “-” in handle to adjust intensity. The function is the same as interface.
9. The system will pause automatically if setting time is over. press any key to continue.
10. Click on “Return” to back to (Fig 1).

### Cautions:

1. Beautician should wear insulated gloves during operation.
2. Do not wear any aglet, both for the beautician and the users.
3. The electrode probe must be fully in contact with skin.
4. Do not make the electrode stay on any part without moving. Keep it moving on skin slowly.
5. Skin around eyes and forehead is very thin, the output intensity should be relatively weakened.
6. Do not contact skin with unipolar when using of 4D circle RV work head otherwise it will produce pulse discharge, gel paste etc. Products like gel maybe evaporated into a mist.
7. Do not let the laser light straightly shine upon the eyes in case of eyes hurt.




## Sprayer

1. Click on “Sprayer” in (Fig 1) to enter interface as (Fig 5).

2. Set work time. About 10 or 20 minutes.

3. Mix 2~3 drip SPA essential oil diluted with purified water into small cup of the oxygen sprayer, water should be about 2 / 3 of the cup, and then tighten the lid tightly. Oxygen sprayer




(Fig 5)

can only with liquid products, other states product will cause pinhole jam.

4. Click on “Start” begin to work. If need to pause during operation, click on “Pause” to stop.

5. Beautician hold the O<sup>2</sup> sprayer, press the control button to release oxygen, and pull backwards to spray water mist. Spray oxygen from top to bottom starting from the forehead to the neck, careful not to spray into the eyes or nose, in the spray oxygen process, the beautician can appease, lift up, stretch wrinkles with the other hand, more helpful for nutrition products penetrate the skin, strengthen skin elasticity, spray rate should depend on the guests' absorb state of oxygen.



6. The system will pause automatically if setting time is over. press any key to continue.

7. Clean the leftovers with hot towel, use warm towel to clean the used working head.

8. Click on “Return” to back to (Fig 1).

### Caution:


1. Avoid spraying into mouth or nose when spraying oxygen.

**NOTICE: If the oxygen sprayer caused blockage or damage while using does not belong to the range of warranty.**

## Cool & Heat

### Cool Therapy

1. Click on “Cool & Heat” in (Fig 1) to enter interface as (Fig 6).
2. The system default is “Cool” therapeutics.
3. Set cool temperature. The system default time is 5°C, click on “-” and “+” of “temperature” to regulate, the range is negative -5°C to 5°C.



(Fig 6)

4. Set work time.
5. Cleaning treated parts, And make zoning, such as: forehead, eyes, cheeks, chin, neck and so on.
6. Do cool treatment according to zone, and adjust appropriate required temperature.
7. Click on “Start” begin to work. If need to pause during operation, click on “Pause” to stop.

### Heat Therapy

1. The system default is “Cool”, Click on “Heat” to select heat therapy.
2. Set working time.
3. Set heat temperature, the system default is 37 °C, click on “-” and “+” to regulate, the range is 37°C to 55°C.
4. Daub oils on treated part, make heat work hand-rod cling to the skin and move slowly.
5. Click on “Start” begin to work. If need to pause during operation, click on “Pause” to stop.
6. The machine will pause automatically if the setting time is over.
7. Click on “Return” to back to (Fig 1).

### Cautions:

1. First use thermal therapy to open pores, select other functions for cleaning or other operations.
2. After the completion of other care, use cool to shrink pores.

## Micro vibration

1. Click on “Double UL” in (Fig 1) to enter interface as (Fig 7).
2. According to treated parts to select big, small ultrasonic work handle. Click on “Large” and “Small” to choose.
3. Select work mode, The system default work mode is P0. click “-” and “+” to choose P1~P5 is pulse wave, P0 is continuous wave.
4. Set intensity. The system default intensity is  $0.5\text{w}/\text{cm}^2$ . click “-” and “+” of intensity to set, the range is  $0.5\text{w}/\text{cm}^2\sim 5.0\text{w}/\text{cm}^2$ .
5. Set work time.
6. Click on “Start” begin to work. If need to pause during operation, click on “Pause” to stop.
7. Apply gel or essence on treatment area, the amount depends on moving flexibility of the ultrasonic head.
8. The machine will pause automatically if the setting time is over.
9. Clean the leftovers with hot towel.
10. Use warm towel to clean and sterilize the used work head.
11. Click on “Return” to back to (Fig 1).

### Cautions:

1. Applied enough gel on skin need treated.
2. Every part treatment should be within 30 minutes.
3. Do not used in the human body heart, throat and eyelids.
4. Do not stay unmoved on any part.
5. The work head can not allowed to be placed under the state of no-load to avoid the damage caused by high heat.



(Fig 7)

---

---

## Touch Screen Calibration

1. Before turn on the system, please click on the screen, don't move the finger, then turn on the Power Switch.
2. The prompts "Err0" is on the screen, "After releasing the key, we enter into the calibration procedure, Operating under the prompts of the screen". Then release your finger now.
3. Operating under the prompts:
  - (1) Press the cross dot in left-top "+".
  - (2) Press the cross dot in right-top "+".
  - (3) Press the cross dot in right-bottom "+".
  - (4) Click on the point (left top).
  - (5) Click on the point (right bottom).
4. Calibration screen is finished, the system enter "input password" interface.
5. Don't input the password, turn off the Power switch, reboot the system.

## Forbidden Group

1. People with heart disease or high blood pressure, or configured cardiac pacemaker.
2. Patients with acute inflammation, asthma, deep vein thrombosis, thyroncus, cancer.
3. People with hemorrhagic disease, trauma or who is bleeding.
4. Pregnant women
5. Medical Plastic parts, or parts with artificial in filling inside.
6. People with metal inside the body(include metallic tooth)
7. Patient with skin inflammation or with edema.
8. People with an abnormal immune system.
9. Numb or insensitive to heat.

## Technical Parameters


Input Voltage: AC220~230V 50/60Hz       110V 50/60Hz

Output Power:  $\leq$ 350VA






Fuse: FL5X20\3.15A with buffer (with spare fuse)

## Packing List

1.Host	1SET	<input type="checkbox"/>
2.BOHR effect C&O handle	1PCS	<input type="checkbox"/>
3.4D-Circle RV handle	1PCS	<input type="checkbox"/>
4.Cool & Heat work handle	1PCS	<input type="checkbox"/>
5.Aqua peeling handle	1PCS	<input type="checkbox"/>
6.Black seal ring(10X1.8)	10PCS	<input type="checkbox"/>
7.Black seal ring(6X1.5)	10PCS	<input type="checkbox"/>
8.Aqua peeling clean cap	10PCS	<input type="checkbox"/>
9.Small ultrasonic work handle	1PCS	<input type="checkbox"/>
10.Big ultrasonic work handle	1PCS	<input type="checkbox"/>
11.Sparyer	1PCS	<input type="checkbox"/>
12.Black seal ring(7.8X1.8)	10PCS	<input type="checkbox"/>
13.User manual	1COPY	<input type="checkbox"/>
14.Power supply line	1PC	<input type="checkbox"/>
15.Liquid bottle(S1-S4)	4PCS	<input type="checkbox"/>
16.Shiny and rehydrate head	1BOX	<input type="checkbox"/>
17.Exfoliation and tightening head	1BOX	<input type="checkbox"/>
18.Left accessories frame	1PCS	<input type="checkbox"/>
19.Right accessories frame	1PCS	<input type="checkbox"/>
20.Handle frame(install right frame)	1PCS	<input type="checkbox"/>
21.Black countersunk hexagonal screw(M4X8)	2PCS	<input type="checkbox"/>
22.Black countersunk hexagonal screw(M4X12)	4PCS	<input type="checkbox"/>
23.Hexagonal screw wrench(M3/M2.5)	2PC	<input type="checkbox"/>
24.Silicone hose with connector	1PCS	<input type="checkbox"/>
25.Silicone massage head	3PCS	<input type="checkbox"/>



## Symbols and Meanings

Warning Symbol	Meaning
 Danger	If used incorrectly will bring the risk of death or serious injury!
 Warning	If used incorrectly will cause possibility ill of death or serious injury!
 Attention	If used incorrectly will cause personnel injury or damage to the goods!
 Danger	Device exists high pressure risk of high voltage output terminals!
 Warning	If used incorrectly will cause personnel injury!

\*Things damage refers to housing, property, livestock and pet damage.

## Contraindications

1. The following symptoms were forbidden to use

- (1) Patients who used by Heart rate regulator, cardiac pacemaker, ECS;
- (2) The use of artificial bone; Patients with silicosis;
- (3) Drunk or high fever patients;
- (4) Patients with higher sensitivity mental disease and epilepsy patients
- (5) Patients with no independent expression consciousness and can not take care of themselves.
- (6) Patients with metal materials in body (including metal teeth), pregnant women, cachexia, malignant tumor.
- (7) Severe hypertension and severe patients with cardiovascular and cerebrovascular disease.
- (8) Patients with renal failure.
- (9) Patients with severe infection



**Danger**

2. The following please use under the guidance of professionals

- (1) Women during menstruation;
- (2) Children;
- (3) Hemorrhagic disease or trauma, inflammation, skin disease, skin infection;
- (4) Patients with feeling numb, no response to cold and heat;
- (5) The immune system has defects, scar formation history, physical disease;
- (6) The people who has taken face-lifting or her body has artificial thing.



**Warning**

## Maintenance

1. Appliances and medical equipment close to the machine may be resulting in operational failure due to interference.
2. Prohibit use of the equipment in the following environment
  - Hot and humid environment
  - Dusty environment
  - Environmental with much smoke or steam
3. If there comes a failure, please shut down the power immediately, and send it to the distributor or manufacturer for maintenance. Do not repair by yourself.

### 4. Exclusion of liability.

**The Company is not responsible for the following circumstances:**

- (1) If fire, earthquake, war occurred;
- (2) The fault of the user, operator error;
- (3) Damage caused by using other method (not included in the use manual).
- (4) Operate without the manual's guide.
5. Do not make the machine close to the wall, keeping a 30CM space around for cooling.

## Safety and Environmental Considerations

### Toxic and hazardous substances reference list Electronic information product pollution control and management measures Table of contents of toxic and hazardous substances or elements

Part Name	Toxic and Hazardous Substance or Element						Remark
	Lead (Pb)	Mercury (Hg)	Cadmium (Cd)	Hexavalent chromium (Cr(VI))	Polybrominated biphenyls (PBB)	Polybrominated diphenyl ethers (PBDE)	
Plastic Board	○	○	○	○	○	○	
Printed Circuit Board	×	○	○	○	○	○	
Power Supply	×	○	○	○	○	○	
Accessory	×	○	○	○	○	○	

#### Declaration:

- The content of the restricted substances in all homogeneous materials of this part is below the limits set by the SJ/T11363-2006 standard.
  - ×
- The content of a homogeneous material in the part of the material is beyond the SJ/T11363-2006 standard requirements. The user's normal use of this product within the time limit of environmental protection, these substances or elements will not occur leakage or sudden change, will not cause the personal, property damage. For such substances or elements, the user shall not dispose of it yourself, please according to government regulations to the relevant departments of the government designated recycling.