


1. Applications

3 Fractional tips with different lengths, more intelligent & accurate for operation

- 1) Ultra Fractional Tip (Short):** Acne, Acne scar, Scar removal, Stretch mark
- 2) Micro-Ablative Tip (Middle):** Wrinkles removal, Pigmentation removal (Freckle, Chloasma, Sun damage)
- 3) Non-Ablative Tip (Long):** Skin Resurfacing



1. Applications

Normal head application:

Surgical cutting (Warts, Nevus, other surgical)



1. Applications

Vaginal head application:

Vaginal tighten, Vulva rejuvenation
Nipples rejuvenation


Vaginal Tighten head



Vulva Treatment




Vulva & Nipples Treatment



The handpiece will be inserted into the vagina area to treat feminine conditions.



When it's inserted, the laser releases thermal stimulation to the ablative area and then to the vagina canal.



After treatment, your vaginal walls starts a natural healing process that regenerates collagen formation in result become firmer & tighter.

2. Advantages

1. Approved by CE certificate
2. 40Watts RF metal tube from America, strong energy and no need maintenance




2. Advantages


3. Why choose RF Metal Tube? Is RF metal tube better than Glass tube?

COMPARISON

RF METAL TUBE



GLASS TUBE



We can provide either RF Metal Tube or Glass Tube (almost US\$2000 Cheaper).

Difference between RF Metal Tube and Glass Tube?


Comparison	Lifetime	Spot size	Downtime	Energy Stability	Maintenance
RF Metal Tube	50,000 hours; 5~8 years	0.12mm	3~4 days	Strong and continuous	Nearly no need maintenance
Glass Tube	20,000 hours; 2 years	0.5mm	7~10 days	Noncontinuous	1~2times/year, high costs

2. Advantages

4. Vaginal rejuvenation system optional,
No consumable, high profits



5. 4 common output patterns and self-designed patterns
by operator, to treat all shapes and areas



6. The most important is that you can not find
on alibaba and the effect is very great!